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G6PDH-GLUCOSE 6 PHOSPHATE DEHYDROGENASE

For *in Vitro* Diagnostic use only**G6PDH QUANTITATIVE KINETIC UV ASSAY on WHOLE BLOOD with MINDRAY BS-300** ∇_{Σ} (20 x 3 mL) 285 tests**REF** NAGP68905

I. INTENDED USE

Glucose-6-phosphate dehydrogenase (G6PDH) is an enzyme in the pentose phosphate pathway. G6PDH converts glucose-6-phosphate into 6-phosphoglucono-d-lactone and is the rate-limiting enzyme of this metabolic pathway that supplies reducing energy to cells by maintaining the level of the coenzyme NADPH. The NADPH in turn maintains the supply of reduced glutathione in the cells that is used to mop up free radicals that cause oxidative damage.

The G6PDH / NADPH pathway is the only source of reduced glutathione in red blood cells (erythrocytes). The role of red cells as oxygen carriers puts them at substantial risk of damage from oxidizing free radicals except for the protective effect of G6PDH/NADPH/glutathione. People with G6PDH deficiency are therefore at risk of hemolytic anemia in states of oxidative stress. Oxidative stress can result from infection and from chemical exposure to medication and certain foods. Broad beans, e.g., fava beans, contain high levels of vicine, divicine, convicine and isouramil, all of which are oxidants.

When all remaining reduced glutathione is consumed, enzymes and other proteins (including hemoglobin) are subsequently damaged by the oxidants, leading to electrolyte imbalance, cross-bonding and protein deposition in the red cell membranes. Damaged red cells are phagocytosed and taken out of circulation in the spleen. The hemoglobin is metabolized to bilirubin (causing jaundice at high concentrations). The red cells rarely disintegrate in the circulation, so hemoglobin is rarely excreted directly by the kidney; but this can occur in severe cases, causing acute renal failure.

Deficiency of G6PDH in the alternative pathway causes the buildup of glucose and thus there is an increase of glycated endproducts. The deficiency also causes a reduction of NADPH which is necessary for the formation of Nitric Oxide (NO).

II. PRINCIPLE

The Glucose-6-Phosphate Dehydrogenase (G6PDH) catalyzes the first step in the pentose phosphate shunt, oxidizing Glucose-6-phosphate (G-6-P) to 6-phosphogluconate (6-PG) and reducing NADP to NADPH.

The increase of absorbance of NADPH, for reduction of NADP, is proportional to the activity of the G-6-PDH in the sample.

III. PRECAUTION FOR USE

1. This product has been formulated for *in vitro* diagnostic use.
2. A proportional variation of the reaction volumes does not change the result.
3. DO NOT mix Reagents from different Production lots.
4. For concentration of G6PDH higher than 3200 U/L, use half sample volume and multiply the result x 2.
5. In addition to the possible risk indications, the Reagent can contain preservatives (as sodium azide or others), which total concentration is lower than the limits mentioned in Dir. 67/548/CEE e 88/379/CEE and following modifications regarding classification, labelling and packaging of dangerous preparations (Reagents). However it is recommended to handle the reagents carefully, avoiding ingestion and contact with eyes, mucous membranes and skin; to use reagents according to good laboratory practice. On the material safety data sheet are detailed the operating procedures for the manipulation of this product. Material safety data sheet should be supplied on request.

ATTENTION!

A) Applications on routine analyzers may be totally different from what developed as manual determination; in addition the procedures are specific for each analyzer.

B) Very deep attention must be given to interfering substances: copper and sulfate are strong inhibitors. Certain drugs and other substances are able to influence levels of G6PDH (see References 2).

C) The reagent must be used ONLY for the intended destinations, by expert and trained people and in according to good laboratory practice.

D) The clinical diagnosis cannot be done correctly using the result of only one test, but have to be done integrating critically the results of different laboratory tests and clinical data.

E)) The reagent contains 6PGD inhibitors, which could form one molecule of NADPH in addition, for one molecule of 6- phosphogluconate formed.

F) Reticulocytes have higher G-6-PDH levels than mature Red Cells; it is not recommended to run the assay after a severe hemolytic crisis, since G-6-PDH may appear falsely elevated.

G) If (□As/min) are very, very low, may be possible increase the sample volume (and of course the Control volume).

H) For the determination of G6PDG on whole blood is necessary to prepare the hemolysate using the kit **REF** NAGB1129 (the procedure is detailed on the relative Information for Use).

I) **G6PDH is very unstable in hemolysates. 20/30 minutes after dilution (see SAMPLE PRETREATMENT NAGREF29) a precipitate may appear, probably due to the biological variability of the patient's sample.**

J) Before to test G6PDH on whole blood it is necessary to determine the following parameter:

- the concentration of Hemoglobin in g/dL, to express G6PDH activity as U/g of Hemoglobin.

K) A lot of factors, as ambient temperature, the working reagent temperature, wash accuracy and the type of spectrophotometer, may affect the tests performances.

L) The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

M) All the precautions normally used in the laboratory must be respected for reagents handling.

All the calibrators and controls must be considered as human sample, so potentially infectious; all the protection actions must be applied to avoid any potential biological risk.

IV. REAGENTS AND MATERIALS PROVIDED

Kit composition:

R1 - BUFFER Good Buffer modified > 20 mmol/L

R2 - NADP NADP >0.19 mmol/L

DANGER

H302-H311-H411-EUH032-EUH208-

P273-P280-P312-P391

REF NAGP668905
1 x 20 mL
20 x 1 mL

R3 - G6P (READY TO USE) G-6-P >0.1 g/L NaN₃ < 0.1% 1 x 40 mL

HAZARD STATEMENTS

- H302** - Harmful if swallowed.
H311 - Toxic in contact with skin.
H411 - Toxic to aquatic life with long lasting effects.
EUH032 - Contact with acids liberates very toxic gas.
EUH208 - Contains Maleimide. It can cause an allergic reaction

PRECAUTIONARY STATEMENTS

- P273** - Do not release into the environment.
P280 - Wear protective gloves / protective clothing / eye protection / face protection.
P312 - Call a POISON CENTER / doctor if you feel unwell
P391 - Collect spilled material.

Contains: sodium azide, Maleimide

MATERIAL REQUIRED BUT NOT PROVIDED

- | | |
|---|---|
| Normal laboratory equipment. | Transparent glass tubes for sample dilution. |
| Micropipettes to deliver from 3 to 1000 µL. | Calibrators and Controls |
| Disposable micropipettes tips. | Spectrophotometer or automatic analyser for Clinical Chemistry. |

AUXILIARY REAGENTS FOR QUALITY CONTROL

To grant the right performances use following kits (see the relative information for use (IFU)):

- | | |
|---|-----------------------|
| - G6PDH-RED CELL LYSING (only for automation) | REF NAGB1129 |
| - TOTAL HEMOGLOBIN | REF NATHEMO395 |
| - G6PDH CONTROLS SET | REF NAG6CON |
| - G6PDH CALIBRATORS | REF NAG6CAL3 |

The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

V. STORAGE AND STABILITY

The Reagents are stable up to the expiry date mentioned on the labels, if closed and stored at 2-8°C in their intact primary container; if not exposed to heat sources and/or pressure variations. In case of damaging of the primary container organize the waste disposal.

VI. PREPARATION OF THE WORKING REAGENT

Dissolve a vial of **R2 - NADP** with 1 mL of **R1 - BUFFER**, mix gently avoid foaming. Let the reagent reach the room temperature before use. Close immediately after handling. The Reagents have to be used correctly, to avoid contamination. An incompetent handling relieves us from any responsibility.

VII. STABILITY OF THE WORKING REAGENT

The working reagent is stable 5 days at 2-8°C.

VIII. SAMPLES

- Whole blood collected with EDTA, heparin or ACD (Acid-Citrate-Dextrose).
- Samples collection in compliance with CLSI (NCCLS) (see References 3).
- Red cell G6PDH is stable in whole blood for 1 week at 2-8°C, but is unstable in Red Cell hemolysate. Freezing of blood is not recommended. (see References 1 et 4).
- Since activity is reported as number of Red Cells or grams of Hemoglobin, these two ones have to be det. before G6PDH assay. The most accurate red cell counts in the time is with erythrocytes in ADC, due to longer integrity of erythrocytes in it than in other ones. (see References 5)

IX. WASTE DISPOSAL

Observe all federal, state and local environmental regulations for waste disposal.

X. ANALYTICAL PROCEDURE FOR WHOLE BLOOD ON MANUAL SPECTROPHOTOMETER

Before to test G6PDH it is necessary to determine:

- 1) the concentration of Hemoglobin in g/dL, to express G6PDH activity as U/g of Hemoglobin.
- Wavelength: 340 nm (334 - 365 nm)
 - Pathlength: 1 cm
 - Reading: against air or distilled water
 - Temperature: 37°C
 - Method: kinetic
 - Reaction: 10 + 2 + 5 minutes
 - Sample/reagent: 1/100/200

Let reagents reach the working temperature before use.

Pipette in a test tube or cuvette so labelled: S: Sample, ST: Calibrators (R1 and R2 **REF NAG6CAL3**)

	S	ST-R1	ST-R2
WORKING REAGENT	1000 µL	1000 µL	1000 µL
Sample	10 µL	---	---
Calibrator R1	---	10 µL	---
Calibrator R2	---	---	10 µL

Mix kindly and incubate for 10 minutes at 37°C. Then add:

R3 - G6P	2000 µL	2000 µL	2000 µL
-----------------	---------	---------	---------

Mix kindly. Exactly 2 minutes after, make the **FIRST** reading of the sample (As1) and calibrators (Ast1-R1) and (Ast1-R2). Repeat the **SECOND** readings after 5 minutes for the sample (As2) and calibrators (Ast2-R1) and (Ast2-R2).

Determine the diff. of absorbance for sample and calibrators:

$$\Delta A_s = A_{s2} - A_{s1} \qquad \Delta A_{st-R1} = A_{st2-R1} - A_{st1-R1} \qquad \Delta A_{st-R2} = A_{st2-R2} - A_{st1-R2}$$

ATTENTION!

The kit is tested for a manual spectrophotometer and for HITACHI, COBAS and MINDRAY systems. The applications on automatic analysers could be completely different by what has been developed as manual determination.

XI. REFERENCE VALUES FOR WHOLE BLOOD (see References 1)

Normal Values G6PDH: 12.1 (7.9 – 16.3) U/g Hb

Values for the newborns may range somewhat higher.

Since the normal values depend on age, sex, diet, geographic area and other factors, each laboratory should establish its own normal values for this procedure.

If, for example, 12.1 U / g Hb corresponds to the average of the reference population and we set it as 100%, the normality values can be expressed in percentage terms as follows: Normal values G6PDH: 65% -135%

XII. CALCULATION

Plot the Δ As values for each Calibrator to give a Calibration Line.

The Calibration Line has to be always repeated for each new lot of Reagent and/or Calibrators.

Plot each Δ As sample value found on the Calibration Line, to obtain its concentration in U/L at 37°C

G6PDH ACTIVITY (U/L - 37°C) Manual Calculation

REM to calculate always ΔA /minute (ΔA /min).

$$\Delta A/\text{min} = (\text{SECOND reading} - \text{FIRST reading}) / 5$$

$$\text{G6PDH (U/L, 37°C)} = \Delta A/\text{min} \times (\text{Total Volume}/\text{Sample Volume}) \times (1/\epsilon \text{ d}) \times 1000$$

$$\text{Total Volume} = (1 + 0,01 + 2) = 3,01\text{mL}$$

$$\text{Sample Volume} = 0,01\text{mL}$$

$$\epsilon = 6,3 = \text{millimolar absorptivity of NADPH at 340nm}$$

$$\text{d} = 1 \text{ cm} = \text{pathlength}$$

$$1000 = \text{Factor to convert activity to Liter}$$

$$\text{G6PDH (U/L, 37°C)} = \Delta A/\text{min} \times (3,01/0,01) \times (1/6,3 \times 1) \times 1000 = \Delta A/\text{min} \times (301 \times 1000) / 6,3$$

$$= \Delta A/\text{min} \times (301000) / 6,3 = \boxed{\Delta A/\text{min} \times 47778}$$

For instance:

$$\text{FIRST reading} = 1,2000$$

$$\text{SECOND reading} = 1,2809$$

$$\Delta A/\text{min} = (1,2809 - 1,2000) / 5 = 0,0809/5 = 0,01618$$

$$\text{G6PDH (U/L, 37°C)} = \boxed{\Delta A/\text{min} \times 47778} = 0,01618 \times 47778 = 773 \text{ U/L (37°C)}$$

G6PDH ACTIVITY (U/gr. Hemoglobin - 37°C) Manual Calculation

Considering the value of Total Hemoglobin (TOTAL Hb) of each sample, expressed as g/dL, apply the formula

$$\text{G6PDH (U/g Hb)} = \frac{\text{G6PDH (U/L, 37°C)}}{\text{TOTAL Hb (g/dL)} \times 10}$$

where "10" is the multiplier that converts g/dL in g/L of Total Hemoglobin (TOTAL Hb).

EXAMPLE

Assay of a sample which had an hemoglobin concentration of 12.47 g/dL and G6PDH (U/L, 37°C) = 773.

$$\text{G6PDH (U/g Hemoglobin)} = \frac{773}{12,47 \times 10} = 6,20$$

XIII. ANALYTICAL PERFORMANCES (Validate on MINDRAY BS300)

The performances of the Reagent **G6PDH - GLUCOSE 6 PHOSPHATE DEHYDROGENASE** have been tested with a MINDRAY BS300 analyzer. The data, while representing the characteristics of the product, could be different for each laboratory and for different analyzers.

Method Limitations: are not know limitations.

Method Linearity: for concentration of G6PDH higher than 3200 U/L, use half sample volume and multiply the result x 2.

Method Sensitivity (LoD): the sensitivity limit, that is the minimum concentration that can be distinguished by zero, is 27 U/L.

Interferences: see References 2.

Interference test criterion: recovery \pm 10% of initial value. No interference found on samples with:

- total bilirubin up to 40 mg/dL;
- lipemia [Intralipid®] up to 4000 mg/dL;
- ascorbic acid up to 50 mg/dL.

Within-run Precision: determined on 20 replications of 2 samples.

The results obtained are following:

Sample	Mean (U/L) \square 2s	CV %
Human 1	191 \square 8	2.1
Human 2	1375 \square 19	0.7

Run-to-run Precision: determined for 5 days with 20 replications for each days, for two samples.

The results obtained are the following:

Sample	Mean (U/L) \square 2s	CV %
Human 1	191 \square 7	1.9
Human 2	1379 \square 19	0.7

Accuracy: a group of 20 sera has been tested using this procedure and using a similar reagent available on the market. The comparison gave these results:

$$\text{Linear regression equation } y = 1.0007x + 15$$

$$\text{Correlation coefficient } r = 0.9997 \text{ n} = 20$$

XIV. REFERENCES (see Italian version)

CONTENT

REF NAGP68905

R1 - BUFFER	1 x 20 mL
R2 - NADP	20 x 1 mL
R3 - G6P	1 x 40 mL
Instruction for use	1 item

In Vitro Diagnostic Medical Device	Temperature limitation	Batch code (AXXX)	Manufacturer
Consult instructions for use	Use by (year/month)	Catalogue number	Do not reuse
Keep dry	Non-sterile	Fragile, handle with care	Keep away from heat