OXALATE

IVD

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BIBLIOGRAFIA

1. Textbook of Clinical Chemistry, Ed. by N.W. Tietz, W.B. Saunders Co., Philadelphia (1999).
2. Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).
3. Tietz of Clinical Chemistry and Molecular Diagnostics, W.B Saunders Co., Philadelphia (2012), (119;154; 581).

CODICE RAMO CND W01010223

OXALATE

IVD

OXALATE QUANTITATIVE COLORIMETRIC TRINDER METHOD ASSAY on URINE
with MINDRAY BS-300 Analyzer

80 tests (1 x 20 mL)

NAOX8850

INTENDED USE

Oxalate in urine is of elevated interest because it could originate from sparingly soluble calcium oxalate; this one could give calcium oxalate crystalluria and the stone formation in the urinary tract, who is considered the most important factor in urolithiasis. Oxalate in urine may increase as an end-product of intermediary metabolism or from dietary sources. A decreased excretion in the urine is associated to hyperglycemia and hyperglycosuria.

An increased excretion could be due to increased ingestion of oxalate rich food or oxalate precursors; an increased one could be also due to metabolic defects or to the absorption of oxalate in several gastrointestinal disorders who give important fat malabsorption.

PRINCIPLE

Oxalate is oxidized to carbon dioxide and hydrogen peroxide by a very specific enzyme, oxalate oxidase.

Hydrogen peroxide reacts in presence of peroxidase (POD) with MBTH (3-methyl-2-benzothiazolinone hydrazone) and DMAB (3-dimethylamino benzoic acid) forming a blue quinone compound. The intensity of colour is proportional to the concentration of OXALATE in the sample and it is read at 590 nm.

Using the standard contained in the kit it is possible to prepare a Calibration Curve to refer.

Plotting on the Calibration Curve absorbance values and concentration for each single sample, may be determined the concentration of each sample.

PRECAUTIONS FOR USE

1. This product has been formulated for in vitro diagnostic use.
2. A proportional variation of the reaction volumes does not change the result.
3. DO NOT mix Reagents from different Production lots.
4. For concentration of Oxalate ≥ 2.00 mmol/L, dilute the sample 1:2 with **distilled water**, repeat the determination and multiply the result by 2.
5. In addition to the possible risk indications, the Reagent can contain preservatives (as sodium azide or others), which total concentration is lower than the limits mentioned in Dir. 67/548/CEE e 88/379/CEE and following modifications regarding classification, labelling and packaging of dangerous preparations (Reagents). However it is recommended to handle the reagents carefully, avoiding ingestion and contact with eyes, mucous membranes and

skin; to use reagents according to good laboratory practice. On the material safety data sheet are detailed the operating procedures for the manipulation of this product. Material safety data sheet should be supplied on request.

ATTENTION!

- A) Applications on routine analyzers may be totally different from what developed as manual determination; in addition the procedures are specific for each analyzer.
- B) Very deep attention must be given to interfering substances: certain drugs and other substances are able to influence levels of Oxalate (see References 2).

C) The reagent must be used ONLY for the intended destinations, by expert and trained people and in according to good laboratory practice.

D) The clinical diagnosis cannot be done correctly using the result of only one test, but have to be done integrating critically the results of different laboratory tests and clinical data.

E) A lot of factors, as ambient temperature, the working reagent temperature, wash accuracy and the type of spectrophotometer, may affect the tests performances.

F) The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

G) All the precautions normally used in the laboratory must be respected for reagents handling.

REAGENTS

Components of the kit:

R1 - BUFFER

(READY TO USE)

Buffer pH 3.1 ± 0.1

> 20 mmol/L

R2 - OXOX, POD

Oxalate Oxidase (Barley)

> 2 KU/L

POD

> 1000 U/L

MBTH

> 0.2 mmol/L

DMAB

> 0.9 mmol/L

Activators, Stabilizers

R3 - PURIFIER

SAMPLE PURIFIER

20 tubes

R4 - DIL

SAMPLE DILUENT CONCENTRATED

1 x 20 mL

R5 - CAL

(READY TO USE)

1 x 5 mL

REF NAOX8850

1 x 20 mL

1 x 2 mL

20 tubes

1 x 20 mL

1 x 5 mL

OXALATE

Oxalic Acid Solution (4.5 mg/dL = 45 mg/L = 0.5 mmol/L)
 NaN₃ < 0.1%

STABILITY: the Reagents are stable up to the expiry date mentioned on the labels, stored at 2-8°C, if closed and kept in their intact primary container; if not exposed to heat sources and/or pressure variations.

In case of damaging of the primary container organize the waste disposal.

AUXILIARY REAGENTS FOR QUALITY CONTROL

To grant the right performances use following kits (see the relative information for use (IFU)):

- CI+OX CALIBRATOR liquid
- OXALATE LOW CONTROL Iyo
- OXALATE HGH CONTROL Iyo
- OXALATE URINE PURIFIER

- REF** OXACAL
- REF** OG6627
- REF** OG6502
- REF** NAOXA75

The calibration curve has to be always repeated at each change of the lot of the Reagent and/or calibrator.

PREPARATION OF THE WORKING STARTER

Dissolve a vial of **R2 - OXOX, POD** with 2 mL of DISTILLED WATER, mix gently until dissolution.

Let the reagent reach the room temperature before use.

Close immediately after handling.

The Reagents have to be used correctly, to avoid contamination.

An incompetent handling relieves us from any responsibility.

STABILITY OF THE WORKING STARTER

The WORKING STARTER is stable 30 days at 2-8°C.

PREPARATION OF THE SAMPLE DILUENT

Add 80 mL of DISTILLED WATER to **R4 - DIL**. Mix gently until complete mixing.

Let the reagent reach the room temperature before use.

Close immediately after handling.

The Reagents have to be used correctly, to avoid contamination.

An incompetent handling relieves us from any responsibility.

STABILITY OF THE SAMPLE DILUENT

The SAMPLE DILUENT is stable 3 months at 2-8°C.

MATERIAL REQUIRED BUT NOT PROVIDED

Normal laboratory equipment.

Micropipettes to deliver from 3 to 1000 µL.

Disposable micropipettes tips .

Transparent glass tubes for sample dilution.

Distilled water, Controls.

Spectrophotometer or automatic analyzer for Clinical Chemistry.

SAMPLES

• No haemolyzed fresh urine, with amount of vitamin C lower than 16 mmol/L, concentration who may affect the test results.

It is recommended that patients refrain from taking excessive amounts of vitamin C or vitamin C-rich-food for at least 48-72 hours prior the urine collection.

Samples collection in compliance with CLSI (NCCLS)

(see References 3): A 24-hour urine sample is collected in a glass or plastic bottle containing 10 mL of conc. Hydrochloric acid. Record the volume in litres.

SAMPLES PRETREATMENT

1. Prepare **SAMPLE DILUENT** as mentioned before.
2. Set up a series of labelled tubes for urine samples and Controls.
3. Pipet 5 mL (or any available volume) of urine samples and Controls into proper labelled tube.
4. Add 2.5 mL of diluted **SAMPLE DILUENT** and 2.5 mL of distilled water into each previous tubes, or any available volume in the same ratio , and mix.
5. Check the pH: it have to be between 5.0 and 7.0. If not, adjust the pH with HCl 1N or NaOH 1N.
6. Set up a series of SAMPLE PURIFIER tubes(**R3 - PURIFIER**) for urine Samples and Controls.
7. Pipet 4 mL each of DILUTED urine Sample and Controls (see point 4. and 5.) to proper labelled SAMPLE PURIFIER tubes (see point 6.) and mix for about 5 minutes by intermittent mixing. We suggest a rotator mixer for mixing.
8. Centrifuge the tubes for 10/15 minutes at 3500 rpm (2600 x g) or filter on laboratory filter paper.
9. Determine the Oxalate concentration picking up the sample in the middle of the supernatant, not from the surface.
10. The Oxalate Standard (**R5 - CAL**) DOES NOT REQUIRE the SAMPLE PREPARATION procedure; it is ready-to-use in the following ANALYTICAL PROCEDURE.

WASTE DISPOSAL

Observe all federal, state and local environmental regulations for waste disposal.

ANALYTICAL PROCEDURE ON MANUAL SPECTROPHOTOMETER

- Wavelength: 590 nm (570-620 nm)
- Pathlength: 1 cm
- Reading: against air or distilled water
- Temperature: 37°C
- Method: end-point
- Reaction: 5 minutes
- Sample/reagent: 1/20/2

Let reagents reach the working temperature before use.

Pipette in 3 test tubes or cuvettes so labelled :

R/B: Blank Reagent, S: Sample and/or Controls, ST: Standard

	R/B	ST	S
R1 - BUFFER	2000 <input type="checkbox"/>	2000 <input type="checkbox"/>	2000 <input type="checkbox"/>
Distilled water	100 <input type="checkbox"/>	----	----
R5 - CAL (do not dilute)	----	100 <input type="checkbox"/>	----
Sample/Control (pretreated)	----	----	100 <input type="checkbox"/>

Mix kindly and incubate at 37°C for 5 minutes. After this one, add:

WORKING STARTER	200 <input type="checkbox"/>	200 <input type="checkbox"/>	200 <input type="checkbox"/>
------------------------	------------------------------	------------------------------	------------------------------

Mix carefully and incubate at 37°C for 5 minutes, waiting the end of the reaction. Read the absorbance of the standard (Ast) and of the sample (As) against the Reagent Blank.

The color is stable for at least 60 minutes.

ATTENTION!

The kit is tested for a manual spectrophotometer and for

OXALATE

IVD

HITACHI, COBAS and MINDRAY systems. The applications on automatic analyzers could be completely different by what has been developed as manual determination.

REFERENCE VALUES (see References 1)

Normal Values Oxalate:

Adults MALES 20 - 60 mg/24 h (0.228 – 0.685 mmol/24 h)
 Adults FEMALES 20 - 55 mg/24h (0.228 – 0.627 mmol/24 h)

Since the normal values depend on age, sex, diet, geographic area and other factors, each laboratory should establish its own normal values for this procedure.

CALCULATION

Use this general formula to calculate the concentration:

Oxalate in mmol/L =

$[As / Astd] \times 0.50$ (Standard conc.) $\times 2$ (dilution factor)]

The Calibration Curve has to be always repeated for each new lot of Reagent.

Oxalate excreted during 24-hour in mmol/24 h =

Oxalate in mmol/L \times Volume of Urine voided in 24-h (in liters)

Oxalate excreted during 24-hour in mg/24 h =

Oxalate in mmol/24 h \times 90 (MW)

ANALYTICAL PERFORMANCES

(validate on MINDRAY BS300)

The performances of the Reagent **OXALATE** have been tested with a MINDRAY BS300 analyzer. The data, while representing the characteristics of the product, could be different for each laboratory and for different analyzers.

Method Limitations: are not know limitations.

Method Linearity: the test is linear up to 2.00 mmol/L.

For concentrations ≥ 2.00 mmol/L, it is recommended to dilute the sample 1:2 with distilled water, test again and multiply the result $\times 2$.

Method Sensitivity (LoD): the sensitivity limit, that is the minimum concentration that can be distinguished by zero, is 0.007 mmol/L.

Interferences: see References 2.

Interference test criterion: recovery $\pm 10\%$ of initial value. No interference found on samples with:

- total bilirubin up to 10 mg/dL;
- haemoglobin up to 600 mg/dL,
- ascorbic acid up to 50 mg/dL.

Within-run Precision: determined on 20 replications of 2 samples.

The results obtained are following:

Sample	Mean (mmol/L) $\pm 2s$	CV %
Human 1	0.233 \pm 0.004	0.7
Human 2	0.943 \pm 0.010	0.5

Run-to-run Precision: determined for 5 days with 20 replications for each days, for two samples.

The results obtained are the following:

Sample	Mean (mmol/L) $\pm 2s$	CV %
Human 1	0.231 \pm 0.005	1.0

Human 2	0.940 \pm 0.013	0.7
---------	-------------------	-----

Accuracy: a group of 20 sera has been tested using this procedure and using a similar reagent available on the market. The comparison gave these results:

Linear regression equation $y = 0.9972x + 0.002$
 Correlation coefficient $r = 0.9997$ $n = 20$

REFERENCES

1. Textbook of Clinical Chemistry, Ed. by N.W. Tietz, W.B. Saunders Co., Philadelphia (1999).
2. Young D.S., Effect of drugs on Clinical Lab. Test, 5th Ed. AACC Press (2000).
3. Titz of Clinical Chemistry and Molecular Diagnostics, W.B Saunders Co., Philadelphia (2012), (119;154; 581).

EDMA (EDMS) CODE 11 02 01 23 00